

**Your partner for advise on and delivery of
Laboratory Animal Equipment!**

ANAEESTHESIA

General Information

The demand for anaesthesia equipment that can be used for small laboratory animals has increased substantially. Following this trend, UNO has been and still is actively involved in the design and manufacturing of complete anaesthesia systems even for the smallest animals. In addition to anaesthesia set-ups with effective evacuation of the anaesthesia gas mixture, we also have monitoring equipment available for Temperature, Capnography, Pulse Oximetry and ECG with Respiratory Monitoring available.

We aim not only to supply the products but make sure that once the products arrive at your facility, the know-how about setting up and using the system is also available! A wide range of products and services is provided to help you to obtain the most suitable set-up.

Pro's and Con's of the three most used type of anaesthesia:

Injectable agents: mostly i.p. or i.m.

- PRO's**
- + Quick and easy to apply (1 injection)
 - + Choice between many products (rat/mouse)
- CON's**
- Due to variation (male/female/breed) or calculated time for procedure, it can be necessary to re-adjust the anaesthesia
 - Injection anaesthesia is IN-FLEXIBLE;
 - extra injection agent if deeper or longer anaesthesia is required.
 - too deep/long anaesthesia when “over injected” or injecting antagonist.

Inhalation anaesthesia: face mask

- PRO's**
- + Flexible anaesthesia “depth” / better control
 - + Also for short(er) procedures
 - + Different gasses (mixtures) possible
- CON's**
- No respiratory control (or pressure control)
 - Adjusting could be relatively slow (if long tubes are used)
 - A gas exhaust system is required.

Inhalation anaesthesia: endotracheal intubation and ventilator

- PRO's**
- + For shorter and longer procedures
 - + Different gasses (mixtures) possible
 - + Full respiratory / -support possible
 - + Quick adjusting when using a **pressure, volume and frequency** driven ventilator.
 - + P.E.E.P. Setting possible (standard AND open thorax)
 - + Economical use of anaesthetic gasses
- CON's**
- Intubation technique is considered to be “difficult” (but with the UNO Intubation Aid it is very easy.)

On the next pages we inform you about the mostly used set-ups available for Inhalation Anaesthesia. But please keep in mind that these are just examples. A large number of adjustments are possible and available.


Just contact us if you need assistance determining what is needed in your facility!

Anaesthesia




General Information


These schemes are just a sample of possible combinations of our anaesthesia products. The dimensions of the separate parts drawn, bear no relation to the real dimensions of the products.



This scheme ABC is an indication of the parts needed to achieve a complete set-up with 1 face mask and 1 gas exhaust unit.



This scheme KMO is an indication of the parts needed to achieve a complete set-up with 1 face mask, 1 induction box and 1 gas exhaust unit.



This scheme TWV is an indication of the parts needed to achieve a complete set-up for controlled anaesthesia with ventilator and intubation aid.



Vaporizers


The Sigma Delta Vaporizer is the latest product from Penlon in a distinguished line of vaporizers of the highest quality and reliability. It delivers accurate concentrations under varying conditions of flow rate and temperature, **particularly at low flows**.



- Selectatec[®], Dräger Plug-In[®], North American Dräger, Cagemount connection
- Superb performance over a wide range of vapour concentration and temperatures, particularly at **low flows**.
- Halothane, Enflurane, Isoflurane, Sevoflurane
- Keyed Filler, Quik Fil[®], or Pour fill
- Efficient Selectatec[®], Dräger, and North American Dräger compatible interlock systems
- Low Body Weight
- **Service Free**; Presumes ten year life requiring no prevetative maintenance service. It is recommended that a service is carried out at ten years. Halothane vaporizers require a five year service.


Technical Specification

| | | |
|-----------------------------|---|-----------------|
| Capacity (ml) | Volume at MAX Mark | 240+/- 10ml |
| | Volume at MIN mark | 35 +/- 10ml |
| | * After draining approx. 60 +/- 10ml of liquid is retained by the wick. | |
| Operating Flow range | 0,2 to 15 L/min | |
| Operating temperature range | 15 to 35°C | |
| Dimensions (wxdxh) in mm | Cagemount | 133 x 158 x 219 |
| | Selectatec Compatible with Interlock | 120 x 190 x 242 |
| | Dräger Plug-In compatible | 100 x 190 x 242 |
| Weight | 5kg | |




The Sigma Delta Vaporizer is available with different connection systems and gasses:


| Vaporizer Type | Isoflurane | | Sevoflurane | | | Enflurane | | Halothane | | |
|--------------------------------------|------------|--------------|-------------|--------------|-----------|-----------|--------------|-----------|--------------|----|
| | 5% | | 8% | | | 5% | 7% | 4% | 5% | 8% |
| Filler Type | Pour | Keyed Filler | Pour | Keyed Filler | Quick Fil | Pour | Keyed Filler | Pour | Keyed Filler | |
| Cagemount Taper mm | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Selectatec Compatible with Interlock | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Dräger plug-In Interlock Compatible | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| North American Dräger | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Selectatec Compatible | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |




Cagemount
Taper



Dräger Plug-In
Interlock Compatible



North American
Dräger



Selectatec Compatible
with Interlock

| Ordering number | Product |
|-----------------|---|
| 18000002 | Isoflurane vaporizer - Cagemount Taper |
| 18000072 | Sevoflurane vaporizer - Cagemount Taper |
| 18000010 | Key- Filler - Isoflurane |
| 18000010-S | Key- Filler - Sevoflurane |
| 18000005 | ISO connector 23mm |

Vaporizers



The flowmeters type SF1, CM2 and SF3 are devices for the supply of medical gases with antistatic and graduated measure tubes and complete with a dosage unit. The body is made of anodized aluminium and the inlet connection at the bottom of the flowmeter is arranged to be connected to existing circuits of gases centralized feeding or to the cylinders complete with pressure regulators. At the outlet of the dosage unit, an anaesthesia vaporizer can be connected or the gas mixture can be supplied directly to the animal through the gas feeding unit. Different constructions can be made on request.



| | SF1 | | CM2 | | SF3 | |
|---|---|---------------|---|----------------|---|----------------|
| Ordering Code | | | | | | |
| Dimensions | | | | | | |
| Height | 300 mm | | 300 mm | | 300 mm | |
| Width | 73 mm | | 113 mm | | 132 mm | |
| Depth | 117 mm | | 117 mm | | 113 mm | |
| Weight | 1,00 Kg | | 1,45 Kg | | 1,9 Kg | |
| Dosage Range | O ₂ | 0,1 - 1 L/min | O ₂ | 0,1 - 2L / min | O ₂ | 0,1 - 1L / min |
| | | | Air | 0,1 - 2L / min | Air | 0,2- 15L / min |
| | | | | | N ₂ O | 0,1 - 1L / min |
| Accuracy | ± 10% read value or ± 0,3 L/min. | | | | | |
| Inlet pressures | 3,5 - 5 bar ± 20% | | | | | |
| Charge loss level, side under pressure (before the adjusting needle valve) | Less than 25ml/min in normal pressure conditions (ISO 5358) | | | | | |
| Charge loss level, low pressure side (after the adjusting needle valves, vaporizer excluded) | Less than 25ml/min. at 30cm H ₂ O. | | | | | |
| Gas Outlet connection | Connical connection* 23 mm F ISO DIN 5356/1 | | --- | | --- | |
| Mixed Gas Outlet connection | --- | | Connical connection* 23 mm F ISO DIN 5356/1 | | Connical connection* 23 mm F ISO DIN 5356/1 | |
| Gasses inlet connection | Tubing connector ø 6 | | Tubing connector ø 7 | | Tubing connector ø 6 | |



*
Connical connection 23mm,
F ISO DIN 5356/1
Ordering code: 180000005



Univentor ; digital controlled anaesthesia unit - combined Flowmeter and Vaporizer

- Designed for small rodents
- Precise control of anaesthetic and air
- Minimised anaesthetic consumption
- Air flow from 50 ml/min up to 999 ml/min
- Connects to mask, anaesthetic chamber or ventilator
- Pre-calibrated for Isoflurane
- Very small foot print and no fixed installation
- User, Animal and Environment friendly
- Audible alarm as end of syringe approaches
- Pusher reverses automatically when syringe is empty
- Easy to use




The Univentor 410 Anaesthesia Unit is designed to control the mixture of liquid anaesthetic and air with the precision required to successfully operate on animals weighing from 20-500 grams.

The Univentor 410 is designed to operate with compressed air **reduced to not more than 1 bar** whereas the Univentor 410-Q has been adapted to work with a pulse-free airpump. Air and anaesthetic is mixed in the vaporizer according the setting and may be delivered into an anaesthetised box, through a mask or to a ventilator.

Even though pre-calibrated for Isoflurane, other anaesthetics may be used taking the various properties into consideration that the unit is calibrated for Isoflurane.

Technical Specification

| | | |
|---------------------------|--|----------------------|
| Dimensions | 120 (W) x 285 (L) x 95 (H)mm | |
| Weight | 1,8 kg | |
| Power supply | 110-240 V, AC 50-60Hz, Battery 12V, 400mA | |
| Drive motor | Pulse free DC Motor with variable speed setting | |
| Fast feed | Pusher movements of 45 mm/min | |
| Pusher Movement tolerance | +/- 0,1mm or +/- 1% of total distance | |
| Max. Pusher Force | 100 N | |
| Concentration tolerance | +/- 0,15% if displayed value | |
| Syringes | 1 glass, gas-tight 10ml syringe with 60mm stroke | |
| Display | 2 x 16 characters | |
| Safety features | Audible alarm and red LED | |
| Min. Liquid Flow rate | 0,4 ml/ hr | |
| Max. Liquid Flow rate | 10 ml/ hr | |
| Min. Air Flow rate | 50 ml/ min | |
| Max. Air Flow rate | 999 ml/ min. | |
| | Univentor 410 - Q | Univentor 410 |
| Min. Air Pressure | 0,3 bar | 0,5 bar |
| Max. Air Pressure | 0,5 bar | 1,0 bar |



When using the Univentor 410-Q you need to have a hydrofobic filter between the pump and the Univentor!




| Ordering number | Product |
|-----------------|---------------------------|
| 180000211 Q | Univentor 410-Q |
| 180000211 | Univentor 410 |
| 180000011-S | Replacement Glass Syringe |
| 180000082 | Pulse free Air Pump |
| 180000080 | Hydrofobic filter |

Univentor



Stands and Trolleys

Our standard **Stand** for mounting the vaporizer and flowmeter is made from black anodised aluminum. Ordering Number: 180000206
Optionally it is available with a holder for the LED holder of the Contrafluran Adsorption filter.



A number of custommade mounts are available as well, like Trolleys or mounts to hang the vaporizer and flowmeter on the wall. Please, contact us if you have specific requirements.



Induction Box

The UNO Induction Box combined with the UNO Gas Exhaust Unit can be used for inducing inhalation anaesthesia with rats and mice. (An anaesthetic gas mixture from a flowmeter/vaporizer should be available.)



The Induction Box is made in 10mm thick red acrylic with an inlet hose connector at bottom level and at the opposite side an outlet $\varnothing 23\text{mm}$ (at the level of the lid).

The special designed lid can be opened by vertical sliding the lid of the Induction Box. When the Induction Box is connected to the UNO Gas Exhaust Unit and the lid is closed, very little air (anaesthetic gas mixture) is being exhausted from the Induction Box because of a relatively high internal resistance. Therefore the animal is optimally exposed to the anaesthetic gas mixture.

As soon as the animal is anaesthetized, the lid of the Induction Box can be slightly slid open. By sliding the lid only a little bit open ($\pm 10\text{mm}$), the internal resistance has gone and the UNO Gas Exhaust Unit is immediately exhausting at maximum capacity, thus preventing the anaesthetic gas mixture to escape from the Induction Box into the room/working area. The Induction Box is thus rapidly emptied from the anaesthetic gas mixture (10 - 15 seconds) and the lid can be taken off to get the anaesthetized animal (mouse or rat) out.

Finally it can be mentioned that the UNO Gas Exhaust Unit can also be connected to an Active Charcoal Filter (for adsorbing the isoflurane) or to an appropriate “in-house” exhaust system.

| Ordering number | Product | Internal Dimensions (L x W x H) | External Dimensions (L x W x H) | Connection Supply hose | Connection Exhaust hose |
|-----------------|--------------------------|---------------------------------|---------------------------------|----------------------------|---------------------------|
| 180000232 | Induction box for mouse | 150 x 900 x 70mm | 170 x 120 x 95mm | $\varnothing 5,4\text{mm}$ | $\varnothing 23\text{mm}$ |
| 180000233 | Induction box for rat | 250 x 130 x 90mm | 270 x 150 x 115mm | $\varnothing 5,4\text{mm}$ | $\varnothing 23\text{mm}$ |
| on request | Induction box for rabbit | 410 x 205 x 268mm | 430 x 225 x 293mm | $\varnothing 5,4\text{mm}$ | $\varnothing 23\text{mm}$ |

Induction
UNO

Open Induction Box



This **OPEN INDUCTION BOX** is meant to be used for situations where a lot of animals have to be treated or handled for a short procedure (blood – or tissue sampling or other short procedures).

The upper part of the OPEN INDUCTION BOX is double walled with holes on the inside and the exhaust tube connections on the outside. The 3 tubes are brought together and all connected to the Exhaust Unit.

When the Exhaust Unit is turned on, an underpressure is created in the double wall. The continuous anaesthetic gas mixture from a Flowmeter / Vaporizer combination enters the lower end of the round Open Induction Box just below the perforated s.s. floor on which the animals are placed.

The gas mixture slowly goes through the perforated floor and stays low because the gas mixture is heavier than air. Of course after building up more volume of the gas mixture, the level rises and finally will be exhausted through the double wall part where the underpressure is.

There is only one important matter: the mouse has to be placed in / taken out **WITH A SLOW MOVEMENT** in order to prevent too much turbulence.

| Ordering number | Product | Internal Dimensions (Ø x H) | External Dimensions (Ø x H) | Connection Supply hose | Connection Exhaust hose |
|-----------------|----------------------|-----------------------------|-----------------------------|------------------------|-------------------------|
| 180000235 | Induction box - Open | 120 x 140mm | 150 x 259mm | ø 23mm | ø 6mm |

Induction



Face Masks

For inhalation anaesthesia without respiration support, UNO has Face masks available;

All models are to be connected to a supply system for an anaesthetic gasmixture (like vaporizer and flow-meter or univentor). The second connection is to a gas exhaust system which safely removes the exhaled/ non used anaesthetic gasses; like the UNO Gas Exhaust unit. The face masks are supplied with extendable exhasut tubes. (40cm - 120cm length).

Face mask for mouse.

- with small nose “cone” specially designed for use with mice
- excellent gas exhaust*
- good cleanability




Face mask for rat.

- with larger nose “cone”
- to be used with rats (and mice depending on type of operation)*.
- gas exhaust and cleanability as for the mouse face mask.

* in combination with the UNO Gas Exhaust Unit.

| Ordering number | Product | Details |
|-----------------|-----------------|--|
| 180000265 | Face mask Mouse | ID Supply: ø 10mm ID Exhaust: ø 20mm Overall dimensions: 60,5 x 28 x 25,25mm Connection Supply hose: ø 4mm Connection Exhaust hose: ø 22mm |
| 180000237 | Face mask Rat | ID Supply: ø 19mm ID Exhaust: ø 29mm Overall dimensions: 81 x 33 x 33mm Connection Supply hose: ø 4mm Connection Exhaust hose: ø 22mm |




Fixator for face masks

The face masks are shaped with a flat bottom but sometimes you want to have it secured to a certain place. For this fixators are available.

- made in Stainless Steel
- to keep the face mask in place, f.i. on a heating plate.

Face mask






Multi Face Mask.

- Allows to have 4 mice kept under anaesthesia simultaneously.
- Each mouse is placed with its snout directly into the face mask*¹.
- The multi face mask is delivered with stoppers that can seal the cones shut.
- On the opposite of the openings for the mouse, the Multi facemask has four (4) connections for face masks with tubing. Thus making it possible to connect 4 Rat face Masks with tubing. It is also possible to have 4 facemasks for Mouse connected if you want to place the animals a bit further apart from each other.

| Ordering number | Product | Details |
|-----------------|---------------------------------------|---|
| 180000136 | Multi face mask | ID Supply: \varnothing 10mm ID Exhaust: \varnothing 20mm Overall dimensions: 264 x 115,5 x 33mm Connection Supply hose: \varnothing 4mm Connection Exhaust hose: \varnothing 22mm |
| 180000136R | Multi face Mask - extended for 4 rats | See Multi face mask plus 4 facemask rat |
| 180000136M | Multi face Mask - extended for 4 mice | See Multi face mask plus 4 facemask mouse |



A face mask for rabbit is also available on request.

UNO GAS EXHAUST UNIT

The UNO Gas Exhaust Unit can be used for the direct-exhaust of waste anaesthetic- and exhaled gas while using Face Masks or Induction Boxes for rat and mouse.

Further the UNO Gas Exhaust Unit can be used for the **indirect-exhaust** of the excessive anaesthetic gas from the circle system of f.i. the UNO Micro Ventilator (UMV) at various volume- and pressure settings of the UNO Micro Ventilator - UMV-03.



In order not to spill the anaesthetic gas mixtures into the working area, and for good functioning of the Gas Exhaust Unit, the supply of the gas mixture to the Face Masks should be set at a maximum of:

- **Face Mask - rat:** - ± 350 ml/min*
- **Face Mask - mouse:** - ± 260 ml/min**
- **UNO Micro Ventilator UMV-03**

The supply of the anaesthetic gas mixture to the UNO Micro Ventilator (UMV) is set at a lower level than the Face Mask, because with the use of the UMV and the correct size intubation tube, there is only little waste of anaesthetic gas mixture***

*** Face Mask Rat:**

- tidal volume: ca. 2,5ml
- breathing frequency during anaesthesia: $\pm 60-70$ /min
- total breathing volume / min: 150 - 175 ml/min
- not all gas supplied is "used" by the rat and as a rule of thumb, the gas supply to the Face Mask is about twice the breathing volume/min, i.e. 300-350ml/min.

**** Face Mask Mouse:**

- tidal volume: ≤ 1 ml
- breathing frequency during anaesthesia: $\pm 110 - 130$ /min
- total breathing volume / min: 110 - 130 ml/min
- not all gas supplied is "used" by the mouse and as a rule of thumb, the gas supply to the Face Mask is about twice the breathing volume/min, i.e. 220 - 260 ml/min.

***** Ventilator UMV-03**

- The supply of the anaesthetic gas mixture to the UMV-03 is set at a lower level than the Face Mask rat or mouse because with the UMV-03 there is only little waste of the anaesthetic gas mixture. It should be slightly higher than the total breathing volume mentioned under Face Mask Rat and Face Mask Mouse plus some extra for maintaining a P.E.E.P. during the ventilation:
 - Rat 250 - 300
 - Mouse 200 - 230

| | |
|-------------------|------------------|
| Ordering number | |
| 180000118 | Gas Exhaust Unit |
| Dimensions | |
| Weight | |
| Power | |
| Inlet Connection | ϕ 23mm |
| Outlet Connection | ϕ 23mm |

Gas Exhaust



Adsorbtion Filters

UNOSORB - FILTER

The UNOSORB Anaesthetic* Gas Filter consists of small grain activated charcoal for optimal adsorbtion of anaesthetic gasses. * Isoflurane, Halothane or Sevoflurane.



- The UNOSORB filter is suitable for Animal Anaesthesia gas removal.
- A disposable canister, containing activated charcoal.
- The UNOSORB canister can be connected directly to the scavenging tubing (f.i. UNO Gas Exhaust Unit).
- The charcoal absorbs organic anaesthetic gasses e.g. Sevoflurane, Isoflurane and Halothane.
- The weight of the canister must be monitored so that it can be replaced before it becomes saturated, once saturated any waste gasses will simply be exhausted into the working environment.




CONTRAFURAN™- FILTER



The patented Contrafluran™ Anaesthetic Gas* Scavenging Filter consists of solid materials, distinguished by their rugged grain structure, extensive surface area and high micro-porosity.

This highly porous internal structure adsorbs efficiently and retains anaesthetic gas components selectively from the exhaled- or un-used anaesthetic gas as it passes through the filter.

The storage capacity of the CONTRAFURAN Anaesthetic Gas Scavenging Filter is ca. 400gr and the flow resistance is with a $\leq 1,5\text{mm wc}$, **VERY LOW!**



The filter can easily be attached to the flowmeter/vaporizer stand or on a separate rail, with the help of the SENSOfluran™, a mount integrated with a visual FILL-LEVEL-CONTROL-UNIT*.

The differently colored LEDs (green, yellow, red) of the Sensofluran Mount indicate the quality of the filtered expired gas and thus the fill level of the filter.




- LED: The filter adsorbs the expired gas and has still sufficient free capacity.
- LED: The capacity of the filter diminishes. A filter change is recommended when the second yellow LED lights up.
- LED: The capacity of the filter is exhausted and should not be used anymore. The used filter must be replaced by a new one.

| Ordering number | Product | Weight | Storage Capacity |
|-----------------|-------------------------------------|-----------------|------------------|
| 180000140 | UNO Sorb Filter | approx. 1.200gr | approx. 200gr |
| 180000138 | Contraflurane Filter | approx. 1.000gr | approx. 400gr |
| 180000139 | LED holder for contraflurane filter | | |

* Due to legal reasons the device must be sent to UNO in a period of 12 months for calibration.

Double 3- and 4-way valves.



To be used for simultaneously changing the direction of both the **SUPPLY** of the anaesthetic gas mixture **AND** the **EXHAUST** of the anaesthetic gas mixture. The risk of making mistakes when the gas supply and gas exhaust have to be changed individually from f.i. the Induction box to the face mask, is reduced substantially by using these valves.



| Ordering number | |
|-----------------|--------------------|
| 180000159 | Double 3-way valve |
| 180000259 | Double 4-way valve |

3-Way Valve:

By turning the knob from the right to the left, both supply and exhaust are being directed from the face-mask to/from the Induction box.




3 way valve




Double 3- and 4-Way Valves


4-Way Valve:



A.. By turning the knob to the left, both supply and exhaust are being directed to/from the induction box.



B.. By turning the knob to the center, both supply and exhaust are being directed to/from the intubation aid.



C.. By turning the knob to the right, both supply and exhaust are being directed to/from the UMV-03.

4 way valve



UNO Micro Ventilator UMV-03

The UNO Micro Ventilator-03 is a mechanical ventilator uniquely designed for use with mice and rats for respiratoric support or totally controlled inhalation anaesthetic.

This revolutionary new concept of mechanical ventilation with an extreme small dead volume, has a proven performance, even with **mice of ≥ 12 gr. bodyweight.**

This unit has even been used for lung transplantation in mice!

The ventilation pattern of the UMV-03 is based on the natural, spontaneous respiratory pattern of rat and mouse, i.e. a **SINUS-form** without a plateau pressure.

The UMV-03 is a **pressure-, volume- and frequency** cycled ventilator which can very easily be integrated in a flowmeter / vaporizer set-up.



Characteristics of the UMV-03:

- ▶ Tidal Volume setting at the ventilator (without counter pressure). 0,1 - 24ml
- ▶ Inspiratory to expiratory ratio 1 : 1 (= sinus form)
- ▶ Effective volume (with counter pressure & endotracheal tube) 0,0ml - 12ml¹
- ▶ Frequency (respiratory rate) 15 - 220/min.
- ▶ Minimum and Maximum pressure setting
- ▶ Pressure settings ≥ 0 mbar - up possible
- ▶ External P.E.E.P.
- ▶ Circle system with CO₂ Absorber
- ▶ Easily to be connected to a flowmeter / vaporizer unit
- ▶ Distance between UMV and animal position upto 1 meter²

| | | | |
|-------------------|----------------------------|-----------------|------------------------------------|
| Ordering number | | Ordering number | |
| 18000023 | UMV-03 UNO Microventilator | 18000009 | Small Artificial Lung |
| Dimensions | | | - used for setting up to UMV-03 |
| Weight | | | |
| Power | | MX00001/1L | CO ₂ Absorber - 1 liter |
| Inlet Connection | ID 5mm | | |
| Outlet Connection | IS 5mm | | |

¹ Depending on the size of endotracheal tube used.

² If distance is more than 1 meter the local situation has to be checked and an extra valve could be necessary.

Intubation Aid for Mice and Rats.

How does it work?

- expose animal to anaesthetics in the induction box for initial anaesthesia
- take animal out of induction box and place “cord” around its front teeth
- slowly pull the animal with its snout/nose into the nose cone with the cord
- secure “cord”
- because the Intubation Aid is connected to the Vaporizer, the animal is **exposed to anaesthetic gas constantly**
- the “**excessive /exhaled gas**” is being disposed off through the connected Gas Exhaust Unit
- after intubation the animal can be connected to the UNO Micro Ventilator.


Benefits:

- *more time for intubation without the need of i.p. anaesthesia; no hurry!*
- *design forces animal in best pose for intubation*
- *lower jaw automatically opens therewith giving clear view on vocal cords*
- *animal can be connected to UMV for controlled anaesthesia immediately after intubation*

Models



- *This model is suitable for mice and small rats ($\pm < 350\text{gr}$)*



- *This model has an exchangeable nose cone, either for mouse or for rats.*
- *Loose nose cones are available. The nose cone can easily be replaced, thus making the intubation aid suitable to be used for mouse and rat.*
- *The design of the aid also makes it usable for bigger rats (extra support of animal.)*

| | |
|-----------------|---|
| Ordering number | |
| 180000014 | Intubation Aid for mice and small rats($\pm \leq 350\text{gr}$) |
| 180000314 | Intubation Aid with Nose cone for Mouse included |
| 180000414 | Intubation Aid with Nose cone for Rat included |
| | |
| 180000324 | Nose cone for Mouse to be used with 180000314 and 180000414 |
| 180000424 | Nose cone for Rat to be used with 180000314 and 180000414 |



Intubation Aid for Mice and Rats.



1.



2.



3.



4.


1. Top part of the intubation stand with tread through “nose-cone”.
2. Mouse placed with snout in nose-cone. Mouse hanging by its front teeth on the tread supported by the back support.
3. Mouse hanging on intubation stand with anaesthetic gas supply- and gas exhaust tubes connected to the intubation stand.
4. Checking the length of the intubation tube (length should be before the bifurcation of the trachea).
- 5. Light source** placed at animal below the vocal cords of the animal. (Do not use a too strong light source because of too much light spreading).
6. Pull the animals tongue a little bit thus enabling you to see the vocal cords (light coming through).
7. Place the intubation tube into the trachea in between the vocal cords.
8. Correct position of the intubation tube with regards to length of the tube and position of bifurcation in trachea.



5.



6.



7.



8.

**The animal is now ready to be connected to the ventilator.
Make sure that the tube maintains its position in the animal!**

Intubation



Controlled Heating System

The UNO Heating devices have been designed specially for the use with small laboratory animals, like mouse and rat during anaesthesia. The UNO Heating systems are powered by a Control Unit.

This Control Unit is available in 3 Versions: **A: Control Unit-01 / CU-01**
B: Control Unit-02 / CU-02
C: Control Unit-MS / CU-MS



Control Unit-MS / CU-MS

Control Unit-02 / CU-02

Control Unit-01 / CU-01

Control Unit-01 and Control Unit-02


The Control Unit 01 has 1, as the Control Unit 02 has 2 connector sockets for one of the below shown heating plates/units and/or a temperature probe. The temperature settings for the Control Unit range from 28°C to 42°C. The connected heating plates/units, heat up till the set temperature on the Control Unit is reached. Once the set temperature is reached the Control Unit maintains this temperature.

The Control Unit 02 does not have an interaction (feedback) between the temperature probe (measured rectal temperature of the animal) and the heating device (set temperature).

The Control Unit-MS does feature this interaction (feedback)!

The **Control Unit-MS** can like the Control Unit-02 be connected to a heating device and a temperature probe. However the interaction (feed back) between the 2 connections allows you **to set the required animals rectal temperature**. The rectal probe measures the animals temperature and feed back is given to the controls of the heating device to keep the animals temperature as constant as possible. A safety feature is built-in to prevent over-heating in case the rectal probe is not in the animal. The heating element will not heat up over 42°C (also an acceptable norm for human skin).


With the Control Unit-MS you have a **self-regulating system** and do not need to manually adjust the temperature settings during the procedure.



Flat Heating plate
 34,5 x 23cm, 12mm thick
 Closed surface for easy cleaning.



Half-Pipe Heating Unit
 for rats during prolonged procedures like cardiac flow measurement.



Tail Heating Unit
 for blood sampling. With or without the half pipe.



Temperature Probe
 for rectal temperature measuring. Mouse or Rat.

| Ordering number | |
|-----------------|---------------------------|
| 180000123 | Control Unit -01 |
| 180000122 | Control Unit -02 |
| 180000124 | Control Unit -MS |
| 180000028 | Flat Plate / 34,5 * 23cm |
| 180000128 | Flat Plate / 18,5 x 11cm |
| 180000027 | Half-Pipe Heating Unit |
| 180000029 | Tail Heating Unit |
| 180000066 | Temperature Probe - Mouse |
| 180000026 | Temperature Probe - Rat |

Warm Water Heating System



The Warm water Heating system contains a warmwaterpump HTP-1500 and a acrylate heating plate. The warm water is set at a certain temperature and is pumped through the canals in the plate. There is no control on the temperature of the animal, the display on the pump provides continous real-time temperature readings. The digital controller uses proprietary software to manage your preset set-pont temperature. Three temperature safety limits guard against possible over-heating.




| | | |
|-----------------|--------------------------------|-----------------------|
| Ordering number | Warm Water Pump HTP-1500 | |
| 180000207 | Water type | Tap water |
| | Reservoir Capacity | 1500 ml |
| | Flow rate (min) | 57 lph |
| | Accuracy watertemperature | +/- 1°C at 42°C |
| | Temperature Set Point Range | 24 to 42°C |
| | Weight | 2,5kg (empty) |
| | Dimensions | 25,5 x 17,8cm |
| | Alarms | Water flow alarm |
| Ordering Number | Heating Plate | |
| 180000021 | Acrylate Heating plate - large | 40 x 30cm |
| Custom made | Acrylate Heating plate - small | 10 x 7,8 20 x 15cm |

Small Animal Retraction System

The UNO retraction systems are a new standard of procedural control. The fully integrated procedural and stabilization system provides precisely controlled retraction at the operative site. The retraction system allows researchers to operate independently, eliminating the requirement for an assistant or improvised assistive devices. The system removes the complexities and distractions of improvised set-ups, allowing the researcher to focus full attention on the surgical procedure.

- Eliminates the need for an assistant
- Reduces setup and procedure time
- Provides superior control
- Improves visualization




The UNO retraction system is unlike any other, providing glove-friendly, operative simplicity. With the touch of a finger, the system can be set, released and manipulated. This is made possible by patented mechanisms in the fixators which use rare earth magnets and push-button spring locks to form a fully integrated small-animal procedure system. All system components are compatible with standard lab instrument cleaning systems.

The System Components

Fixators- The Heart of the System; Fixators hold retractors. They contain rare earth magnets that attach anywhere on the animal table. A push-button top operates multi-level locking jaws that grip wires or elastomers firmly. They can be adjusted linearly or rotationally with finger-tip release. Complete depression of the push-button allows for insertion of the wires or elastomer while partial depression of the push-button allows for adjustment. Fixators can also lock onto other lab components that require stabilization in the surgical setting, such as neurological or rectal probes or anaesthesia delivery systems.

The magnetic field is well controlled within the fixator to minimize interference with sensitive instrumentation.

In fact, at a distance of just a few centimeters the effects of these small but powerful magnets are completely gone.



Retraction




Small Animal Retraction System


UNO Retraction System - Mouse

Contains:

- 1pcs, 20cm x 30cm table with window
- 6pcs, 3cm fixators
- 10pcs, 10cm wire retractor handles
- 10cs, 2 meter roll elastomer
- 10pcs, assorted retractor tipe, two of each style
- User documentation (english)



Magnetic Fixator with Spring Lock
3cm tall, Two locking jaws
ACD-001




Mouse Table 20 x 30cm
With 13 x 10cm Window
ACD-003

| | |
|--|---|
| | Retractor tip Sharp - 1mm (10 per pack) ACD-009 |
| | Retractor tip Blunt - 1mm (10 per pack) ACD-010 |
| | Retractor tip Blunt - 2,5mm (10 per pack) ACD-011 |
| | Retractor tip Blunt - 5mm (10 per pack) ACD-012 |
| | Retractor tip Blunt - 7,5mm (10 per pack) ACD-013 |
| | Retractor Wire- 10cm (10 per pack) ACD-005 |
| | Retractor Wire- 14cm (10 per pack) ACD-006 |
| | System Elastomer (2 meter roll) ACD-007 |
| | Retractor Tip- Sample Package (2 each style, 5 sizes per pack) ACD-008 |


UNO Retraction System - Rat

Contains:

- 1pcs, 25cm x 35cm table with window
- 6pcs, 6cm fixators
- 10pcs, 14cm wire retractor handles
- 10cs, 2 meter roll elastomer
- 10pcs, assorted retractor tipe, two of each style
- User documentation (english)



Magnetic Fixator with Spring Lock 6cm
tall, Four locking jaws
ACD-002



Rat Table 25 x 35cm
With 25 x 20cm Window
ACD-004

Retraction




Small Animal Retraction System

Animal Tables- Animal Tables are made of ferro-magnetic stainless steel and incorporate a window for compatibility with body temperature maintenance systems. Table are available in a variety of sizes and shapes to suit a range of animals.



Retractors- The retractor wires are formed from light, flexible stainless steel. This simple approach produces an economical instrument that performs its job with minimal clutter, while providing superb control and feel. Retractor wires can either be used on their own or retract or manipulate tissue or retractor tips can be fitted to provide a wide range of retractor styles. Retractor wires engage and lock into the fixator jaws at the desired elevation allowing both linear and rotational adjustment over the entire length of the wire. One end of the wire is formed in a loop to provide grip for a gloved hand, and the business end is formed into a right angle hook that can retract on its own or lock to a range of retractor tips.




The wires can also be reformed to create restraints or engage other lab apparatus that need to be secured in the surgical field. Wires are available in 10 and 14cm lengths for mouse and rat procedures. Retractor tip widths range from 1mm to 7.5mm, plus a single point, sharp tissue hook. Each tip can be locked onto a wire for static retraction, or affixed to an elastomer to form a dynamic retractor.

| Ordering number | |
|-----------------|-------------------------|
| 18000055 | Retraction System RAT |
| 18000056 | Retraction System MOUSE |

Retraction



The new **UNO Cold Light Source** is designed to meet the current requirements for a cold light source during surgical procedures.



Advantages

- New power - LED technology
- Colour temperature 6300°K (Pure white)
- No cooling ventilator required
- Narrow beam lens: 13° View Angle
- With rechargeable battery with up to 10 hours battery life at maximal light strength.
- Power consumption: 700mA, Power max: 4 Watt
- Luminous flux max. 100lm - dimmable
- Power requirement 110 - 240V, 50/60Hz
- Very long service life of Power LED'S
- Very small footprint (14 x 12 x 5cm (W x D x H))

In order to maintain a distance between the CLS and the illuminated object, the CLS is equipped with two (2) flexible 60cm long arms (custom lengths available on request).

| | |
|-----------------|-------------------------|
| Ordering number | |
| 180000901 | Cold Light Source - CLS |



Replacement parts

| | Ordering Number | Description |
|---|------------------|---|
|  | 180000024 | Silicone tubing Ø 4 x 7mm - sold per meter |
|  | 180000025 | Silicone tubing Ø 6 x 12mm - sold per meter |
|  | 180000040 | Red PVC Tubing Ø 6 x 8,6mm - sold per meter |
|  | 180000234 | Corrugated extensible Tube ID 22mm, - sold per 40 -120cm |
|  | 180000034 | Ribbed (Harmonica) Tubing - sold per 40cm |
|  | 180000045 | Connector 22M / 22F - ID 22mm- sold per piece |
|  | 180000046 | Connector 22M - 22M - sold per piece |
|  | 180000052 | Tube valve - sold per piece |
|  | 180000047 | 3 way valve (PP/PE) - sold per piece |
|  | 180000044 | Parallel Y-Connector - sold per piece |

Pulse Oximeter VET Handheld Pulse Oximeter

The 2500A Veterinary Oximeter is a universal tool with unmatched versatility and cost-effectiveness for all veterinary monitoring and research applications. Compact and easy to use, this Pulse Oximeter is proven accurate for pulse rates up to 450 beats per minute (bpm) and is ideal for monitoring during surgeries.

The 2500A Pulse Oximeter has an audible pulse indicator and large LED displays that are visible in low-light situations. The tricolor perfusion indicator provides immediate feedback to assess pulse quality. This information is useful to determine if repositioning of the sensor is necessary and requires little training to interpret.

Sensors

Purelight Vet sensors produce the pure light spectrum which eliminates variations in readings from animal-to-animal and sensor-to-sensor.

2000SL - Clip-on sensor for tongue applications, paw pads, and well-vascularized areas.

2000T - For placement on the underside base of the tail or other well-vascularized surfaces.

2000SA - Wrap sensor for placement on the toe (large animal) or base of the tail or foot (small animal).



Key features

- * Durable - extremely rugged construction
- * Extensive Memory - 72 hours data storage
- * Easy to Use - Simple two button operation
- * Compact Size - 7 x 13,8 x 3,2cm
- * Flexible Operation - operates on AA batteries and rechargeable batteries

Pulse Oximeter Specifications

| | |
|---------------------------------|---|
| Oxygen Saturation Display Range | : (%SpO ₂) 0 -100% |
| Pulse Range Display Range | : 18 - 450 beats per minute (BPM) |
| Accuracy | : Blood oxygen Saturation (%SpO ₂ ± 1 SD) 70 - 100% ± 2 digits. Pulse Rate: 18 - 450 bpm ±3% ± 1 digit |
| Measurement Wavelengths | : Red 660 Nanometers Infrared 925 Nanometers |
| Dimensions | : 7 x 13,8 x 3,2cm (w x h x d) |
| Weight | : 213 gr with alkaline batteries, 233 grams with NiMH batteries |
| Temperature | : Operating -20°C to +50°C Storage -30°C to +50°C |
| Humidity | : Operating 10 - 95% noncondensing Storage 10 - 95% noncondensing |
| Power requirements | : Four 1,5V AA-size alkaline batteries |

Pulse Oximeter

SENSOR PLACEMENT

The **2000SL LINGUAL CLIP SENSOR** is convenient for spotchecks and for monitoring during recovery while the animal is immobile. However the sensor is held in place by light spring pressure that could allow it to become dislodged with movement. Further, over time even light spring compression may interfere with blood flow resulting in signal loss and need to reposition the sensor. The 2000SL is most easily applied to the rear foot, but the front foot and tail can also be used. For mice, the sensor can be best placed high into the groin.



The **2000T TRANSFLECTANCE SENSOR** also suitable for continuous monitoring, the 2000T transreflectance sensor is the smallest probe and can be applied to the rear of front feet, or to the ventral surface of the tail of the rat. Secured with adhesive tape.



Alternatively the sensor can be attached to the surface of an operating or imaging table and the foot simply taped into position over it.

The **2000SA SMALL ANIMAL WRAP SENSOR** can be easily secured and less prone to accidental displacement, making it an ideal option for continuous monitoring during long surgical or other procedures. Adhesive tape, cohesive bandage or a Posey Wrap can be used to assure that the sensor diodes are aligned directly opposite each other and that excessive pressure is avoided. The 2000SA should not be applied to hair-covered or highly pigmented areas.



| Ordering number | |
|-----------------|-------------------------------------|
| 180000143 | Pulse Oximeter 2500 VET incl 2000SL |
| 2000SL | Sensor Lingual clip Sensor |
| 2000T | Transflectance sensor |
| 2000SA | Small Animal Wrap |

* Above mentioned sensor placement recommendation are taken from the article Pulse Oximetry for Rodents by Dr. George A. Vogler, DVM. A copy is available on request.

* More detailed information on sensor placement and reported possible sensor sites are available on request.

Capnograph AMP for Rats and large Animals (not suitable for use with Mice!)

CAPNOTRUE® AMP MAINSTREAM CO₂ /SPO₂ MONITORS

Advanced and reliable capnograph and pulse oximeter combined in a single monitor.



Key features of CapnoTrue® AMP

- Mainstream CO₂ measurement with the IRMA™ CO₂ analyzer
- Warm-up time: < 10s full specification
- Direct measurement without time delay
- Small, light-weight and shock-resistant: the IRMA™ CO₂ analyzer weighs less than 30 g
- Adult/paediatric, infant/neonatal, and a **for rat modified** IRMA™ airway adapter
- Plug and measure technology
- IRMA™ airway adapters with non-condensing light transmission XTP™ window
- No occlusion by water or mucus possible
- Maintenance and calibration-free technology
- Full range of key technology accessories
- Wide range of high-quality SpO₂ sensors

Delivery Package of CapnoTrue® AMP

- Mainstream device
- IRMA™ CO₂ analyzer
- Reusable SpO₂ sensor
- CapnoTrue® power supply (EU and UK plug)
- Silicone protective cover
- USB data cable
- IRMA™ airway adapter (adult/paediatric)
- Instruction manual
- 1 Li-ion rechargeable battery
- PC software
- 4 batteries (AA)

The CapnoTrue AMP Mainstream Monitor can be used in combination with the UNO MICROVENTILATOR - UMV-03.

The CapnoTrue AMP and ASP - CO₂/SpO₂ monitors including the IRMA™ CO₂ analyzer, ISA™ CO₂ analyzer and SoftCap Sensors are classified and certified as class IIb products.

Capnograph



Capnograph AMP for Rats and large Animals (not suitable for use with Mice!)

TECHNICAL DATA

| Specification | | |
|---------------------------------|--|--|
| Measurement range | EtCO ₂ and FICO ₂ | 0-15% |
| | SpO ₂ | 0-100% |
| | Respiration rate | 0-150 breaths/min |
| | Pulse rate | 20-300 beats/min |
| Accuracy | EtCO ₂ and FICO ₂ | +/- (0.2 vol % + 2% of reading) +/- (0.3 vol % + 4% of reading) incl. interfering gases |
| | SpO ₂ | +/- 2% (70 - 100%) |
| | Respiration rate | +/-1 digit at 60 breaths/min |
| | Pulse rate | +/-1 digit (up to 100/min) or +/- 1% (> 100/min) |
| Display | | |
| Characteristics | Active OLED colour graphic display, 262,000 colours, 240 x 320 dots | |
| Displayed data | End-tidal CO ₂ and inspired CO ₂ in vol %, kPa or mmHg, oxygen saturation, respiration rate, pulse rate, capnogram, plethysmogram and short term trend | |
| Indicators | Signal strength and signal quality, pulse amplitude, battery status, alarm mute, pulse tone mute, neonatal mode, time | |
| Trend information | | |
| Long term trend | up to 400 hours | |
| Short term trend | 15 min/1h/6h | |
| Language versions | | |
| | English | (additional language versions on request) |
| Environmental conditions | | |
| Operation | 0 - 50°C, 15 - 95% r.H. (non-condensing), 60 - 120 kPa (excl. Li-ion battery) | |
| Storage | -30 - 70°C, 10 - 95% r.H. (non-condensing), 60 - 120 kPa (excl. Li-ion battery) | |
| Classification | | |
| Product class | IIb (in accordance with MDD 93 / 42 / EEC) | |
| Safety | Class of protection II / type BF | |
| Construction | IPX1 | |
| Standards | EN 60601-1:2006, EN 60601-1-1:2001, EN 60601-1-2:2007, EN 60601-1-8:2004, DIN EN 60529, EN ISO 21647:2004, EN ISO 9919:2005 | |
| Miscellaneous | | |
| Dimensions | (LxWxD) 150 mm x 75 mm x 35 mm | |
| Weight | < 400 g (complete device with batteries) | |
| Power supply | 90-260 VAC/6 VDC, 50-60 Hz, 4 AA batteries, 2.5 Ah Li-ion battery | |
| Battery capacity | up to 6 hours continuous operation | |
| Communication interface | USB 2.0 | |

CapnoTrue MG for Rats and large Animals (not suitable for use with Mice!)

CAPNOTRUE® MG-AA/CO₂/SPO₂ MONITORS

High-performance and versatile anaesthetic agent monitoring.



Key features of CapnoTrue® MG

- Innovative micro-optic technology
- Direct mainstream measurement without time delay
- Compact, robust and ultra-light multigas analyzer
- Warm-up time <20 seconds full specification
- Maintenance and calibration free technology
- Self-explanatory, ergonomic operating function facilitate intuitive operation
- The colour information display, as well as the simple information structure, support quick decisions and a rapid user reaction in critical situations
- Leading-edge power management with standard alkaline batteries or Li-Poly batteries or medical power supply (or combined)
- Two years warranty.

Delivery Package of CapnoTrue® MG

- CapnoTrue MG Multigas/SpO₂ Monitor
- IRMA™ AX+ Analyzer
- Instruction manual
- 1 Li-ion rechargeable battery
- PC software
- 4 batteries (AA)
- IRMA™ Airway Adapter
- Reusable SpO₂ sensor
- Silicone protective cover
- USB data cable
- Manual

The CapnoTrue MG Monitor can be used in combination with the UNO MICROVENTILATOR - UMV-03.

The CapnoTrue MG monitors including the IRMA™ CO₂ analyzer, ISA™ CO₂ analyzer and SoftCap Sensors are classified and certified as class IIb products.

Capnograph MG

TECHNICAL DATA

Display

| | | |
|----------------------|--|--|
| Parameters displayed | Numerical | End-tidal (et) CO ₂ -, N ₂ O- and agent concentrations, inspired (Fi) CO ₂ -, N ₂ O- and agent concentrations, oxygen saturation (SpO ₂), Respiration Rate (RR), Pulse Rate (PR) |
| | Graphical | Capnogram and trends of numerical data (15 min/1 h/6 h) |
| | Indicators | Signal Strength and signal quality, pulse amplitude, battery status, alarm mute, pulse tone mute, storage status, real-time mode, neonatal mode, time. |
| Characteristics | Active OLED colour graphic display, 262000 colours, 240 x 320pixel (42mm x 56mm) | |

Capnography and anaesthetic agent measurement specifications

| | | |
|-------------------|---|---|
| Measurement range | etCO ₂ and FICO ₂ | 0-15% |
| | FIN ₂ O | 0-100% |
| | Hal, Iso, Enf | 0-8% |
| | Sev | 0-10% |
| | Des | 0-22% |
| | Respiration rate | 0-150 l/min |
| Accuracy | EtCO ₂ and FICO ₂ | +/- (0.2 vol % + 2% of reading) +/- (0.3 vol % + 4% of reading) incl. interfering gasses |
| | N ₂ O | +/- (2 vol % + 2% of reading) +/- (3 vol % + 5% of reading) incl. interfering gasses |
| | Hal/Iso/Enf/Sev/Des | +/- (0,15 vol % + 5% of reading) +/- (0,2 vol % + 10% of reading) incl. interfering gasses |
| | Respiration rate | +/-1 digit |
| Warm-up time | < 20 seconds full specification | |

Pulse Oximetry specifications

| | | |
|-------------------|------------------|---|
| Measurement Range | SpO ₂ | 1 -100% |
| | Pulse Rate | 20 -300 1/min |
| Accuracy | SpO ₂ | +/- 2% (70 bis 100%) |
| | Pulse Rate | +/- 1 digit (up to 100 1/min) or +/- 1% (> 100 1/min) |

Trend Information

| | |
|------------------|-----------------|
| Long term trend | up to 150 hours |
| Short term trend | 15 min/1h/6h |

Alarms

| | |
|--------|---|
| Limits | Adjustable limits for all numerical parameters except for MAC |
| Alerts | Audible and visual alarms (complies with EN60601-1-8) |

Storing Data

| | |
|-------------------------|---|
| Communication interface | USB 2.0 |
| Data memory capacity | up to 150 hours |
| Real-time mode | Visualisation and storage of numerical parameters on a computer every 4 seconds |
| Computer Software | CapnoTrue® MG PC Software for data download and real-time mode |

Power supply

| | |
|-----------------|--|
| Battery | Working time with full functionality approx. 4,5 hours, 4 alkaline batteries (AA/LR6/AM3/MN1500/Mignon), 1,5V |
| Li-Poly battery | Working time with full functionality approx. 7 hours, Li-Poly battery, Model No CT-2500, 3.7V, 2500 mAh, charging time approx 5 hours |
| AC power supply | Model No. FW 7660M/06, medical power supply with option for country-specific input plug, input: 100-240 V AC/50-06Hz / 250 mA, output: 6 V DC/1,4A |

CapnoTrue MG for Rats and large Animals (not suitable for use with Mice!)

TECHNICAL DATA

| Specification | |
|---------------------------------|--|
| Environmental conditions | |
| Operation | 10 - 40°C, 15 - 95% r.H. (non-condensing), 60 - 120 kPa (excl. Li-poly battery) |
| Storage | -20 - 70°C, 10 - 95% r.H. (non-condensing), 60 - 120 kPa (excl. Li-poly battery) |
| Classification | |
| Product class | I Ib (in accordance with MDD 93 / 42 / EEC) |
| Safety | Class of protection II / type BF - Type and degree of protection against shock |
| Construction | IPX1 (with silicone protective cover) |
| Standards | EN 60601-1:2006, IEC 60601-1-1:2001, IEC 60601-1-4-2000; IEC 60601-1-8:2006; ISO 21647:2004; ISO 9919:2005; ISO5356-1:2004; DIN EN 1789:2007, EN846:1996 |
| Miscellaneous | |
| Dimensions | 150 mm x 75 mm x 35 mm |
| Weight | < 400 g (complete device with batteries) |


With its ultra-compact, light and easy to handle design, the CapnoTrue®MG is the perfect flexible and mobile monitor for identifying and quantifying the five most important anaesthetic agents as well as other gases and parameters:

- Halothane, Isoflurane, Enflurane, Sevoflurane and Desflurane
- N₂O
- etCO₂, FICO₂
- Oxygen saturation
- Respiration and Pulse Rate

By direct measurement in the mainstream, there is no time delay in the measuring data.

Reliable automatic agent identification and quantification

The IRMA AX+ Analyzer is equipped with state-of-the-art NDIR technology with up to 9-channel gas type analysis in the 4-10µm range and offers reliable agent identification and quantification even in gas mixtures. It weighs less than 25 g.






Probes to be used with rat

Probes to be used with larger animals

Special designed IRMA Airway Adaptor with a very small dead space.

UMV-03 - Ventilator




The Phasein CO₂ measuring unit is placed over the special IRMA Airway Adaptor. The mainstream measurements are taken through the small round window in the adaptor.

| Ordering number | |
|-----------------|--------------------------------------|
| 180000170 | Capnotrue AMP |
| 180000171 | Capnotrue MG (Multigas) |
| 180000172 | Adjusted Probe to be used with Rats |
| 180000173 | Probe to be used with larger Animals |



Surgical Monitoring Mouse and Rat Made Easy



Heart Rate * SpO2 * Temperature * ECG * Warming

Better Results

Get better study results while improving surgery quality & survival rate.

Superior Data

Ultra low-noise, high-resolution ECG, SpO2 & Respiration.

Easy to Use

Durable stainless steel surgical platform AND intuitive touch display.

Stable Prep

Maintain body temp. & monitor vital signs during surgery.

Physiological Parameters

ECG

- Ultra low-noise, high resolution
- Simultaneous Lead I, II & III
- 24-bit sigma-delta AD converters

Heart Rate

- Real-time numeric display
- Up to 999BPM
- Acquired from ECG waveform

Pulse Oximetry

- Ultra low-noise, high resolution
- Simultaneous SpO2 display
- 18-bit sigma-delta AD converters

Respiration

- Ultra low-noise, high resolution
- Based on thoracic impedance
- Signal acquired through ECG electrodes

Breath Rate

- Real-time numeric display
- Up to 300BrPM
- Acquired from respiration waveform

Core Temperature

- 0,1°C monitoring accuracy
- Mouse-specific rectal probe

Warming

- Electronic closed-loop control
- Intelligent zone heating
- Platform temperature control with 0,1°C resolution


Intuitive Touch Interface

25,5cm Touch screen

Capacitive touch screen works with surgical gloves

Ergonomic Display

Waveforms and numeric data presented in an easily readable format



Comments

Apply notes, tags, and observations to mark events

Screenshot

Take quick snapshots to capture interesting data

Record and Export

Supports many analysis packages through CSV export

Heated Surgical Platform

Accessories

Frame designed for magnetic accessories

Expansion Modules

Four expansion ports support analog output and SpO2 modules

Durable



Easy to clean stainless steel work surface and electrodes

Warming Zone

Homothermic heater maintains core temperature.

ECG Electrodes

Surface-mounted mouse & rat electrodes enable easy operation with low noise



MouseMonitor™ also suitable for RAT

The MouseMonitor™ is a compact tabletop vital signs monitor that displays ECG and respiration waveforms as well as heart rate, breath rate and core temperature. The unique integrated pad incorporates ultra-low noise, high resolution ECG electronics and a homothermic heating pad with a durable surgical steel operating surface that supports magnetic accessories. The included 10.1” touchscreen display ergonomically presents waveform and numeric data in an intuitive and legible layout, enabling you to quickly record and export your acquired data. The MouseMonitor™ is compatible with most DAQ systems when connected to the optional analog output module.

ECG


- Uses non-invasive pad mounted or external ECG electrodes
- Simultaneous Lead I, Lead II and Lead III ECG
- Works in supine or prone position, head in either direction
- Electrically isolated ECG bioamplifier
- High-resolution 24 bit sigma-delta AD converters

Respiration

- Uses non-invasive pad mounted ECG electrodes

Temperature Control and Monitoring

- Electronic closed-loop control
- Efficient, electrically heated warming zones
- Pad temperature control with resolution of 0.1°C
- Compact tabletop design without circulating water



Specifications

Heated Surgical Platform

| | |
|----------------|--|
| ECG Electrodes | 4 mouse limb electrodes 4 rat limb electrodes External electrodes with 1,5mm DIN jacks |
| Heater | Electronic heating with multiple zone controle |
| Temperature | Adjustable, 25-42 °C |
| Audible Alarm | Yes |
| Size | 25,4 x 30,4cm |
| Weight | +/- 3 kg |

Display Unit

| | |
|------------|---|
| Display | 25cm color, LED backlit, capacitive touch |
| Resolution | 1280 x 800 pixels |
| Storage | 8GB |
| Power | 100-240V AC adapter |

Rectal Temperature Probe

| | |
|-------------|--------------------------------|
| Probe Type | Mouse rat or neonatal specific |
| Sensor Type | Copper-Constantan thermocouple |
| Accuracy | +/- 0,1°C after calibration |

Measurements

Numeric

- Heart Rate 60-999BPM
- Respiration Rate 15-400BrPM
- Core Temperature 25-50°C
- Pad Temperature 25-42°C

Waveform


- ECG Lead I
- ECG Lead II
- ECG Lead III
- Respiration

| Ordering number | |
|-----------------|--|
| 180000300 | Mousemonitor |
| 180000301 | Analog Output Mode |
| 180000302 | Pulse Oximetry/SPO2 Module |
| 180000303 | Platinum External Needle Electrodes |
| 180000304 | Stainless Steel External Needle Electrodes |

UNO Euthanasia Unit

During euthanasia of rodents using the UNO Euthanasia Unit, the following benefits are achieved:

- The animals can stay in their cage, because the UNO Euthanasia Unit includes a transparent polycarbonate “container”, in which the complete rat- or mousecage can be placed.
- First, Carbogen gas is introduced at the bottom of the cage, by means of a PLC controlled valve. Carbogen is a gas mixture of 95% O₂ + 5% CO₂. This causes a high concentration O₂ in the cage and in the container, while the 5%-concentration CO₂ does not cause any reaction at the CO₂ receptors in the animal lungs. The pO₂ level in the blood becomes maximal.
- After 1 min. the Carbogen flow is stopped and 100% CO₂ is introduced into the cage. In this way the high concentration O₂ is slowly replaced by CO₂. This relatively slow replacement of O₂ by CO₂ is responsible for the reduction of stress for the animal.




If mice and rats are being exposed directly to 100% CO₂, the stress level, measured according to EEG, EKG and proportionally, is very high. The article “Carbon Dioxide euthanasia in rats;” handles about the 3 above mentioned criteria in detail. This article can be downloaded on our website:

<http://www.unobv.com/afbeeldingen/Carbon%20Dioxide%20in%20rats.pdf>

The UNO Euthanasia Unit includes:

- A polycarbonate box with lid, in which the mice and rat cages can be placed.
Dimensions: 50 x 75 x 25cm (lxbxh) (other dimensions custom made available!)
- Tube for the gassupply. This tube can be put between the bars of the wirelids into the cage.
- A valve in the macrolon box to connect the unit to an exhaust system.
- The s.s. trolley has a tray on which the gas cylinders can be placed.
- A switchboard for the supply of the gas (mixture) O₂ and CO₂, according a PLC pre-set program.




One of the advantages of the mentioned PLC pre-set programs is that the procedures for euthanasia in rats and mice can be easily described in a protocol for the users of the unit. The UNO Euthanasia Unit contains a program that uses 2 different gasses in order to, during euthanasing mice and rats, reduce the animal's stress to a minimum.

| | |
|-----------------|--------------------------|
| Ordering number | |
| 180000100 | Euthanasia Unit complete |

Euthanasia




Different configurations and set-ups for euthanasia units are available on special request. Below are some examples.



A Euthanasia Unit with two separate euthanasia boxes that can each hold a cage. This unit was equipped with two separate control with each just one program. Also a safety to prevent unauthorized changes made in the flowmeter settings was integrated.



A Euthanasia Unit to be placed on a table without the trolley is also possible.



Euthanasia cabinets in different configurations. For ordering information, please contact us.

Oxygen Concentrator 5.0

- Modular components set up
- Easily accessible filter and fuses
- USB – interface
- Modern design and simple operation
- New, innovative measuring sensor
- Integrated flow adjustment
- 30.000 hours warranty for all functional parts
- Whisper silent
- Microprocessor-controlled



Accessories

- Instructions for use
- Angular connector for humidifier
- Humidifier, refillable
- Nasal cannula 2m and 5m

| Technical Data | |
|-------------------------------|---|
| MDD classification | IIa |
| Operating Power | 230 V 50Hz |
| Operating Temperature | operation +5° to +40°C storage -20° to +70°C |
| Sound level | < 40 dB (A) |
| Power consumption | 295 W |
| Dust and Fine Filter | in the rear of the device |
| Fuses | mains T 3.15 A H 250V internal T 1.0 A L 250V |
| Weight | 15kgs |
| Dimensions (HxLxW) | 600 x 290 x 400mm |
| Warranty | 30.000 operating hours (max. 5 years) according to terms of warranty |
| Flow adjustment | flowmeter 0 - 5 lpm |
| O ₂ -concentration | 1 to 3 lpm 95% -3% 3 to 4 lpm 93% +/- 3% 4 to 5 lpm 85% +/- 3% |

| | |
|-----------------|-------------------------|
| Ordering number | |
| 180000099 | Oxygen Concentrator 5.0 |

Recovery / Nursery / Warming

For improving the temperature of weak animals or animals recovering after a surgical procedures, we have a flexible recovery system available. The system is based on the cage type IVS and can be used as a compact system with one or two cages, placed on desk, table of other available space or delivered in a moveable rack.



A silicone rubber mat is placed on the heating plate at the bottom of the cage. This mat can be changed easily in order to minimize contact with an other mouse or rat in the cage after a period of recovery/nursery.

The heating plate on the bottom of the cage can be taken out through a slot at the short side of the cage. This slot is closed and the closure keeps the heating plate in place. The heating system is controlled by a heating Control Unit which is positioned on the top of the rack. The temperature range is 28°C- 42°C and can be easily adjusted.



We have two options for the cover of the cage.

1. A perforated polycarbonate plate covering the total opening of the cage. In this perforated plate is a small lid that can be opened to be able to put the animal (or water and feed) into the cage. With this system the warmth will stay longer in the cage because the total openings are smaller.
2. A conventional wire lid. With this wire lid the animal is still in contact with the heating system via the silicone rubber mat. The warmth does not stay in the cage too long; it will go out through the wire lid.



| | |
|--------------------|--|
| Ordering number | |
| Please contact us. | |

Recovery



UNO Roestvaststaal BV

PO Box 15 - NL 6900AA Zevenaar

T: +31 316 524451

E: info@unobv.com - [Http: www.unobv.com](http://www.unobv.com)